

Product datasheet for **RC223329L3V**

NABC1 (BCAS1) (NM_003657) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NABC1 (BCAS1) (NM_003657) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NABC1
Synonyms:	AIBC1; NABC1; PMES-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_003657
ORF Size:	1752 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223329).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003657.1
RefSeq Size:	3475 bp
RefSeq ORF:	1755 bp
Locus ID:	8537
UniProt ID:	O75363
Cytogenetics:	20q13.2
MW:	61.7 kDa



[View online »](#)

Gene Summary:

This gene resides in a region at 20q13 which is amplified in a variety of tumor types and associated with more aggressive tumor phenotypes. Among the genes identified from this region, it was found to be highly expressed in three amplified breast cancer cell lines and in one breast tumor without amplification at 20q13.2. However, this gene is not in the common region of maximal amplification and its expression was not detected in the breast cancer cell line MCF7, in which this region is highly amplified. Although not consistently expressed, this gene is a candidate oncogene. [provided by RefSeq, Apr 2016]