

Product datasheet for **RC223344L3V**

RGAG4 (RTL5) (NM_001024455) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	RTL5
Synonyms:	6430402L03Rik; MAR5; MART5; RGAG4; SIRH8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PSI00092)
Tag:	Myc-DDK
ACCN:	NM_001024455
ORF Size:	1707 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC223344).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001024455.1
RefSeq Size:	4356 bp
RefSeq ORF:	1710 bp
Locus ID:	340526
UniProt ID:	Q5HYW3
Cytogenetics:	Xq13.1



MW: 64.5 kDa