

Product datasheet for **RC223346L3V**

PRPF38B (NM_018061) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PRPF38B (NM_018061) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PRPF38B
Synonyms:	NET1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018061
ORF Size:	1638 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223346).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018061.1 , NP_060531.1
RefSeq Size:	3239 bp
RefSeq ORF:	1641 bp
Locus ID:	55119
UniProt ID:	Q5VTL8
Cytogenetics:	1p13.3
Domains:	PRP38
Protein Pathways:	Spliceosome



[View online »](#)

MW: 64.3 kDa

Gene Summary: May be required for pre-mRNA splicing.[UniProtKB/Swiss-Prot Function]