

Product datasheet for **RC223364L4V**

ZNF618 (NM_133374) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	ZNF618
Synonyms:	FP13169; NEDD10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_133374
ORF Size:	2583 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC223364).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_133374.2 , NP_588615.2
RefSeq Size:	9133 bp
RefSeq ORF:	2586 bp
Locus ID:	114991
UniProt ID:	Q5T7W0
Cytogenetics:	9q32



MW: 95.4 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]