

Product datasheet for RC223366L3

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

PDLIM5 (NM_001011516) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: PDLIM5 (NM_001011516) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: PDLIM5

Synonyms: ENH; ENH1; L9; LIM

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clor

Sequence:

The ORF insert of this clone is exactly the same as(RC223366).

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_001011516

ORF Size: 642 bp





PDLIM5 (NM_001011516) Human Tagged Lenti ORF Clone - RC223366L3

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001011516.2</u>, <u>NP 001011516.1</u>

RefSeq Size: 1971 bp
RefSeq ORF: 645 bp
Locus ID: 10611
UniProt ID: Q96HC4
Cytogenetics: 4q22.3

Protein Families: Druggable Genome

MW: 23.8 kDa

Gene Summary: This gene encodes a member of a family of proteins that possess a 100-amino acid PDZ

domain at the N terminus and one to three LIM domains at the C-terminus. This family member functions as a scaffold protein that tethers protein kinases to the Z-disk in striated muscles. It is thought to function in cardiomyocyte expansion and in restraining postsynaptic growth of excitatory synapses. Alternative splicing of this gene results in multiple transcript

variants. [provided by RefSeq, Jan 2012]