

## Product datasheet for **RC223377L3V**

### **MURF3 (TRIM54) (NM\_187841) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	MURF3 (TRIM54) (NM_187841) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MURF3
Synonyms:	MURF; MURF-3; muRF3; RNF30
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_187841
ORF Size:	1074 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223377).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_187841.1</a>
RefSeq Size:	1644 bp
RefSeq ORF:	1077 bp
Locus ID:	57159
UniProt ID:	<a href="#">Q9BYV2</a>
Cytogenetics:	2p23.3
MW:	40.1 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene contains a RING finger motif and is highly similar to the ring finger proteins RNF28/MURF1 and RNF29/MURF2. In vitro studies demonstrated that this protein, RNF28, and RNF29 form heterodimers, which may be important for the regulation of titin kinase and microtubule-dependent signal pathways in striated muscles. Alternatively spliced transcript variants encoding distinct isoforms have been reported. [provided by RefSeq, Jul 2008]