

## Product datasheet for **RC223410L2V**

### **GABA A Receptor alpha 6 (GABRA6) (NM\_000811) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | GABA A Receptor alpha 6 (GABRA6) (NM_000811) Human Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | GABRA6   |
| Mammalian Cell Selection: | None   |
| Vector:                   | pLenti-C-mGFP (PS100071)   |
| Tag:                      | mGFP   |
| ACCN:                     | NM_000811  |
| ORF Size:                 | 1359 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC223410).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_000811.1</a>  |
| RefSeq Size:              | 1732 bp  |
| RefSeq ORF:               | 1362 bp  |
| Locus ID:                 | 2559   |
| UniProt ID:               | <a href="#">Q16445</a>   |
| Cytogenetics:             | 5q34   |
| Domains:                  | Neur_chan_memb, Neur_chan_LBD  |
| Protein Families:         | Druggable Genome, Ion Channels: Cys-loop Receptors, Transmembrane  |



[View online »](#)

|                          |  |
|--------------------------|--|
| <b>Protein Pathways:</b> | Neuroactive ligand-receptor interaction  |
| <b>MW:</b>               | 51.02 kDa  |
| <b>Gene Summary:</b>     | GABA is the major inhibitory neurotransmitter in the mammalian brain where it acts at GABA-A receptors, which are ligand-gated chloride channels. Chloride conductance of these channels can be modulated by agents such as benzodiazepines that bind to the GABA-A receptor. At least 16 distinct subunits of GABA-A receptors have been identified. [provided by RefSeq, Jul 2008] |