

Product datasheet for **RC223410L4V**

GABA A Receptor alpha 6 (GABRA6) (NM_000811) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GABA A Receptor alpha 6 (GABRA6) (NM_000811) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GABA A Receptor alpha 6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_000811
ORF Size:	1359 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223410).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000811.1
RefSeq Size:	1732 bp
RefSeq ORF:	1362 bp
Locus ID:	2559
UniProt ID:	Q16445
Cytogenetics:	5q34
Domains:	Neur_chan_memb, Neur_chan_LBD
Protein Families:	Druggable Genome, Ion Channels: Cys-loop Receptors, Transmembrane



[View online »](#)

Protein Pathways:	Neuroactive ligand-receptor interaction
MW:	51.02 kDa
Gene Summary:	GABA is the major inhibitory neurotransmitter in the mammalian brain where it acts at GABA-A receptors, which are ligand-gated chloride channels. Chloride conductance of these channels can be modulated by agents such as benzodiazepines that bind to the GABA-A receptor. At least 16 distinct subunits of GABA-A receptors have been identified. [provided by RefSeq, Jul 2008]