

## Product datasheet for **RC223412L3V**

### FOXC2 (NM\_005251) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	FOXC2 (NM_005251) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FOXC2
Synonyms:	FKHL14; LD; MFH-1; MFH1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_005251
ORF Size:	1503 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223412).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005251.1</a>
RefSeq Size:	1506 bp
RefSeq ORF:	1506 bp
Locus ID:	2303
UniProt ID:	<a href="#">Q99958</a>
Cytogenetics:	16q24.1
Domains:	FH
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

**MW:** 53.5 kDa

**Gene Summary:** This gene belongs to the forkhead family of transcription factors which is characterized by a distinct DNA-binding forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in the development of mesenchymal tissues.  
[provided by RefSeq, Jul 2008]