

Product datasheet for **RC223427A1V**

Human Heterogeneous Nuclear Ribonucleoprotein (A1 like) (HNRNP A1L2) (NM_001011725) AAV Particle

Product data:

Product Type:	AAV Particles
Tag:	Myc-DDK
Symbol:	Heterogeneous Nuclear Ribonucleoprotein (A1 like)
Mammalian Cell	None
Selection:	
Vector:	pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC223427 representing NM_001011725
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGTCTAAGTCAGCGTCTCCAAAAGAGCCGAACAGCTGAGGAAGCTCTTCATTGGAGGGTTGAGCTTTG
AAACAACCTGATGAGAGCCTGAGGAGCCATTTTGAGCAATGGGGAACGCTCACAGACTGTGTGGTAATGAG
AGATCCAAACACCAAGCGCTCCAGGGGCTTTGGGTTTGTACATATGCCACTGTGGAGGAGGTGGATGCA
GCTATGAATACAACGCCACACAAGGTGGATGGAAGAGTTGTGGAACCAAGAGAGCTGTCTCCAGAGAAG
ATTCTCAAAGACCAGGTGCCCACTTAAGTGTGAAAAAGATATTTGTTGGTGGCATTAAAGAAGACTGA
AGAACATCACCTAAGAGATTATTTGAACAGTATGGAATAATTGAAGTAATTGAAATCATGACTGACCGA
GGCAGTGGCAAGAAAAGGGGCTTTGCCTTTGTAACTTTGACGACCATGACTCCGTGGATAAGATTGTCA
TTCAGAAATACCTACTGTGAAGGGCCACAACCTGTGAAGTTAGAAAAGCCCTGCCAAAGCAAGAGATGGC
TAGTGCTTCATCCAGCCAAAGAGGTCGAAGGGGTTCTGGAACTTTGGTGGTGGTCTGGAGATGGTTTC
GGTGGGAATGACAACCTTTGGTCTGGAGGAACTTCAGTGGTCTGGTGGCTTTGGTGGCAGCTGTGGT
GTGGTGGATATGGTGGCAGTGGGGATGGCTAATGGATTTGGTAATGATGGAAGCAATTTGGAGGTGG
TGGAAGCTACAATGATTTTGGCAATTACAACAATCAGTCTTCAAATTTGGACCATGAAGGAGGAAAT
TTTGGAGGCAGAAGCTCTGGCCCTATGGCGGTGGAGGCCAATACTTTGCAAAACCACAAAACCAAGGTG
GCTATGGCGTTTCCAGCAGCAGCAGTAGCTATGGCAGTGGCAGAAGATT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:	>RC223427 representing NM_001011725 Red=Cloning site Green=Tags(s) MSKSASPKEPEQLRKLFIGGLSFETTESLRSHFEQWGLTDCVVMRDPNTRSRGFGFVTYATVEEVDA AMNTTPHKVDGRVVEPKRAVSREDSQRPGAHLTVKKIFVGGIKEDTEHHLRDYEYQYKIEVIEIMTDR GSGKKRGFAFVTFDDHDSVDKIVIQYHTVKGHNCEVRKALPKQEMASASSQRGRGSGNFGGGRGDGF GGNDNFGRGGNFSGRGGFGGSCGGGGYGGSGDGYNGFGNDGSNFGGGGSYNDFGNYNQSSNFGPMKGGN FGGRSSGPYGGGGQYFAKPQNQGGYGVSSSSSYGSGRRF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001011725
ORF Size:	960 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001011725.1</u>
RefSeq Size:	2224 bp
RefSeq ORF:	963 bp
Locus ID:	144983
UniProt ID:	<u>Q32P51</u>
Cytogenetics:	13q14.3
MW:	34 kDa