

Product datasheet for **RC223428**

PTPN20B (PTPN20) (NM_001042365) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PTPN20B (PTPN20) (NM_001042365) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: PTPN20
Synonyms: bA42B19.1; bA142I17.1; CT126; PTPN20A; PTPN20B
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC223428 representing NM_001042365
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTGGACAGCCAGAGGCCCTTCAGAAGACAGGTGGAGCAGTGAGGATGAGGAGGCTGCAGGGCCAT
CACAGGCTCTCTCCCTCTACTTTCTGTTCAACATCATGGATATAGTGGCCAAATGAGAGAACAACGTT
CTGGCATGGTTCAAACGAAGGAGCAGTATCACTTTTGTACGATATTGTGCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC223428 representing NM_001042365
Red=Cloning site Green=Tags(s)

MWTARGPFRRDRWSEDEEAAGPSQALSPLL SVQH HGYSGPNERTTFWHGSNEGAVSLLLR YCA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001042365
 ORF Size: 192 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001042365.3
RefSeq Size:	1885 bp
RefSeq ORF:	195 bp
Locus ID:	26095
Cytogenetics:	10q11.22
Protein Families:	Druggable Genome, Phosphatase
MW:	7 kDa
Gene Summary:	The product of this gene belongs to the family of classical tyrosine-specific protein tyrosine phosphatases. Many protein tyrosine phosphatases have been shown to regulate fundamental cellular processes. The encoded protein appears to be targeted to sites of actin polymerization. A pseudogene of this gene has been defined on chromosome 10. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2014]