

Product datasheet for **RC223445L4V**

RNF35 (TRIM40) (NM_138700) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	RNF35
Synonyms:	RNF35
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138700
ORF Size:	687 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC223445).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_138700.2
RefSeq Size:	997 bp
RefSeq ORF:	690 bp
Locus ID:	135644
UniProt ID:	Q6P9F5
Cytogenetics:	6p22.1



MW: 25.9 kDa

Gene Summary: This gene encodes a member of the tripartite motif (TRIM) protein family. The encoded protein may play a role as a negative regulator against inflammation and carcinogenesis in the gastrointestinal tract. Alternatively spliced transcript variants that encode different protein isoforms have been described. [provided by RefSeq, Feb 2014]