

Product datasheet for **RC223456L1V**

Adropin (ENHO) (NM_198573) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Adropin (ENHO) (NM_198573) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Adropin
Synonyms:	C9orf165; UNQ470
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_198573
ORF Size:	228 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223456).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_198573.2
RefSeq Size:	1093 bp
RefSeq ORF:	231 bp
Locus ID:	375704
UniProt ID:	Q6UWT2
Cytogenetics:	9p13.3
Protein Families:	Transmembrane
MW:	7.7 kDa



[View online »](#)

Gene Summary: Involved in the regulation of glucose homeostasis and lipid metabolism.[UniProtKB/Swiss-Prot Function]