

## Product datasheet for **RC223484L3V**

### LSM11 (NM\_173491) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	LSM11 (NM_173491) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LSM11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_173491
ORF Size:	1080 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223484).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_173491.2</a> , <a href="#">NP_775762.1</a>
RefSeq Size:	2563 bp
RefSeq ORF:	1083 bp
Locus ID:	134353
UniProt ID:	<a href="#">P83369</a>
Cytogenetics:	5q33.3
MW:	39.3 kDa


[View online »](#)

**Gene Summary:**

Component of the U7 snRNP complex that is involved in the histone 3'-end pre-mRNA processing (By similarity). Increases U7 snRNA levels but not histone 3'-end pre-mRNA processing activity, when overexpressed. Required for cell cycle progression from G1 to S phases. Binds specifically to the Sm-binding site of U7 snRNA.[UniProtKB/Swiss-Prot Function]