

Product datasheet for **RC223498A1V**

Human Speedy protein C (SPDYC) (NM_001008778) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Speedy protein C (SPDYC) (NM_001008778) AAV Particle
Tag:	Myc-DDK
Symbol:	Speedy protein C
Synonyms:	Ringo2; RINGOC
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC223498 representing NM_001008778 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCTCTGGGCCATTCTGAGCTCGGGTCACCCTGCCCCATCTCCATCTCCTATGAGATGAGTGACTCCC
AAGACCCACCACTTCCCCTGTAGTTACCAACCAGGTAGAGCTGGGGGGCTGCAGCCGGCAAGGTGGGGG
CAATGGGTTCCCTCCGTTTTCCGCAGCACCAGGAGGTCCAGGCCTTCTCAGCCTTCTGGAGGACAGTTTT
GTCCAGGAATTCCTCTCCAAAGACCCCTGCTTCCAGATTCAGATAAGTATCTCCTGGCCATGGTGCTGG
TCTACTTCCAGCGCGCCACCTGAAGCTCAGCGAGTATACCCACAGCAGCCTGTTCTTGGCCCTGTACCT
TGCAAACGACATGGAGGAGGACCTGGAGGGCCCCAAATGTGAGATTTTTCCATGGGCCCTGGGAAAAGAT
TGGTGTTTACGAGTGGGAAAATTCCTGCACCAGAGGGATAAGCTTTGGGCACGGATGGGTTTCCGGGCTG
TTGTGAGCCGCCAGTGCTGTGAGGAGGTGATGGCAAAGGAGCCATTCCACTGGGCTTGGACTCGGGACCG
GCGCCCCACCATGGTGGGTTCCAGAGGTCTGTCCACAGGTCCCTGTTGCCTTCCCCGGGGCCCTGGC
CTCTCGCCGCCCACTGTTCCCCTGTGGTTTGGCCAGCACTGCAGCAGCCACTTGCTTAAGCCTGTGT
CATCAAGTGCCCTTCTCTGACCTCTGAATGTCATCGCCCTCCCTCCCAAATTATCTCTCAAGGTCAA
AAACGCCTGGGGTGGGACTTTCTCATCGTCTTGCTTCCCGAGATGCAACTGGAACAGGCACCTACTCC
CTCCGCATCTTCCAAAGCCTCCGGCACGCCCTGGGCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC223498 representing NM_001008778 Red=Cloning site Green=Tags(s) MLWAIPELGSPCPIISISYEMSDSQDPTTSPVVTQVELGGCSRQGGNGFLRFRQHQEVQAFLSLLEDSF VQEFLSKDPCFQISDKYLLAMVLVYFQRAHLKLSEYTHSSLFLALYLANDMEEDLEGPKCEIFPWALGKD WCLRVGKFLHQDKLWARMGFRAVVSRCCEEVMAKEPFHWAWTRDRRPHHGGVQVCPQVPVRLPRGPG LSPPHCSPCGLPQHCSHLLKPVSSKCPSLTSECHRPPSQNYLSRVKNAWGGDFLIVLPPQMQLPEPGTYS LRIFPKPPARPGH TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001008778
ORF Size:	879 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001008778.1
RefSeq Size:	882 bp
RefSeq ORF:	882 bp
Locus ID:	387778
UniProt ID:	Q5MJ68
Cytogenetics:	11q13.1
MW:	33 kDa