

## Product datasheet for **RC223504L4V**

### Methyltransferase like 6 (METTL6) (NM\_152396) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Methyltransferase like 6 (METTL6) (NM_152396) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Methyltransferase like 6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_152396
ORF Size:	852 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223504).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_152396.2</a>
RefSeq Size:	2499 bp
RefSeq ORF:	855 bp
Locus ID:	131965
UniProt ID:	<a href="#">Q8TCB7</a>
Cytogenetics:	3p25.1
Protein Families:	Druggable Genome
Protein Pathways:	Androgen and estrogen metabolism, Histidine metabolism, Selenoamino acid metabolism, Tyrosine metabolism



[View online »](#)

**MW:** 33.1 kDa

**Gene Summary:** Probable methyltransferase.[UniProtKB/Swiss-Prot Function]