

Product datasheet for **RC223505A1V**

Human C17orf98 (NM_001080465) AAV Particle

Product data:

Product Type: AAV Particles

Tag: Myc-DDK

Symbol: C17orf98

Mammalian Cell None

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC223505 representing NM_001080465
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGTACCTGAGCGAGTGTGCGCTGCGACTGGAGAAAGGCTTTATCTTGGACGGGTGGCTGTGAGCA
CCGCTGCCCCGCGCTTATGGGCGCTCTAGGCCAAGCTGTGGTCGGCGATTCCGCCCTACAACGCGCAGCA
GGACTACCACGCCCGCAGCTACTTCCAGAGTCACGTGGTTCGCCCCCTTTGCGGAAAACCTGATCAGGAT
CATGGCGGTACAGGAAGGGATGGCTGGATAGTAGACTACATCCACATCTTCGGGCAAGGACAAAGATACC
TCAACAGGAGAAATTGGGCAGGACAGGGCATTCTCTCCAGCAGGTGACTGGGCATGACCACTACAATGC
TGATCTGAAACCGATCGATGGGTTCAATGGAAGGTTTGGCTACCGCAGGAACACTCCAGCCCTCCGTCAG
AGCACGTCCGTCTTCGGAGAGGTCAACCACTTCCCTCTGTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC223505 representing NM_001080465
Red=Cloning site Green=Tags(s)

MAYLSECLRLLEKGFILDGVAVSTAARAYGRSRPKLWSAIPPYNAQQDYHARSYFQSHVVPPLLRKTDQD
HGGTGRDGWIVDYIHFQGGQRYLNRRNWAGTGHSLQQVTGHDHYNADLKPIDGFNGRFGYRRNTPALRQ
STSVFGEVTHFPLF

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2

ACCN: NM_001080465



ORF Size:	462 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001080465.1
RefSeq Size:	465 bp
RefSeq ORF:	465 bp
Locus ID:	388381
UniProt ID:	A8MV24
Cytogenetics:	17q12
MW:	17.4 kDa