

Product datasheet for **RC223511A1V**

Human **RAB27A** (NM_183236) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human RAB27A (NM_183236) AAV Particle
Tag: Myc-DDK
Symbol: RAB27A
Synonyms: GS2; HsT18676; RAB27; RAM
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC223511 representing NM_183236
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCTCTGATGGAGATTATGATTACCTCATCAAGTTTTAGCTTTGGGAGACTCTGGTGTAGGGAAGACCA
GTGTACTTTACCAATATACAGATGGTAAATTTAACTCCAAATTTATCACACAGTGGGCATTGATTTAG
GGAAAAAGAGTGGTGTACAGAGCCAGTGGGCCGGATGGAGCCACTGGCAGAGGCCAGAGAATCCACCTG
CAGTTATGGGACACAGCAGGGCAGGAGAGGTTTCGTAGCTTAACGACAGCGTTCTTCAGAGATGCTATGG
GTTTTCTTCTACTTTTTGATCTGACAAATGAGCAAAGTTTCCTCAATGTCAGAACTGGATAAGCCAGCT
ACAGATGCATGCATATTGTGAAAACCCAGATATAGTGCTGTGTGGAAACAAGAGTGATCTGGAGGACCAG
AGAGTAGTGAAAGAGGAGGAAGCCATAGCACTCGCAGAGAAATATGGAATCCCCTACTTTGAAACTAGTG
CTGCCAATGGGACAAACATAAGCCAAGCAATTGAGATGCTTCTGGACCTGATAATGAAGCGAATGGAACG
GTGTGTGGACAAGTCTGGATTCTGAAGGAGTGGTGCATCAATGGTCATGCCTCTACGGATCAGTTA
AGTGAAGAAAAGGAGAAAGGGGCATGTGGCTGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC223511 representing NM_183236
Red=Cloning site Green=Tags(s)

MSDGDYDYLKFLALGDSGVGKTSVLQYTDGKFNSKFITTVGIDFREKRVVYRASGPDGATGRGQRIHL
QLWDTAGQERFSLTTAFFRDAMGFLLLFDLTNEQSFLNVRNWISQLQMHAYCENPDIVLCGNKSDLEDQ
RVVKEEEAIALAEKYGIPYFETSAANGTNIQAIEMLLDLIMKRMERCVDKSWIPEGVVRNNGHASTDQL
SEEKEKGACGC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_183236
ORF Size:	663 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_183236.1
RefSeq Size:	2772 bp
RefSeq ORF:	666 bp
Locus ID:	5873
UniProt ID:	P51159
Cytogenetics:	15q21.3
MW:	24.7 kDa