

Product datasheet for **RC223513**

SPFH2 (ERLIN2) (NM_001003791) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: SPFH2 (ERLIN2) (NM_001003791) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: SPFH2
Synonyms: C8orf2; Erlin-2; NET32; SPFH2; SPG18
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC223513 representing NM_001003791
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTCAGTTGGGAGCAGTTGTGGCTGTGGCTTCCAGTTTCTTTGTGCATCTCTTCTCAGCTGTGC
ACAAGATAGAAGAGGGACATATTGGGGTATATTACAGAGCGGTGCCCTGCTGACTTCGACCAGCGGCC
TGGTTTCCATCTCATGCTCCCTTTCATCACATCATATAAGTCTGTGCAGACCACACTCCAGACAGATGAG
GTGAAGAATGTACCTTGTGGACTAGTGGTGGTGTGATGATCTACTTTGACAGAATTGAAGTGGTGAAC
TCCTGGTCCCGAACGCAGTGTATGATATAGTGAAGAACTATACTGCTGACTATGACAAGGCCCTCATCTT
CAACAAGATCCACCACGAACTGAACCAGTTCTGCAGTGTGCACACGCTTCAAGAGGTCTACATTGAGCTG
TTTGACTGGAATGATTTTTCCAGGAATCTTCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC223513 representing NM_001003791
Red=Cloning site Green=Tags(s)

MAQLGAVVAVASSFFCASLFSAVHKIEEGHIGVYYRGGALLTSTSGPGFHLMLPFITSYKSVQTTLQTDE
VKNVPCGTSGGVMIFYDRIEVVNFVLPNAVYDIVKNYADYDKALIFNKIHHELNQFCSVHTLQEVYIEL
FGLENDFSQESS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8040_c08.zip



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001003791

ORF Size: 456 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001003791.3](#)

RefSeq Size: 1538 bp

RefSeq ORF: 459 bp

Locus ID: 11160


UniProt ID: [O94905](#)

Cytogenetics: 8p11.23

MW: 16.8 kDa

Gene Summary: This gene encodes a member of the SPFH domain-containing family of lipid raft-associated proteins. The encoded protein is localized to lipid rafts of the endoplasmic reticulum and plays a critical role in inositol 1,4,5-trisphosphate (IP3) signaling by mediating ER-associated degradation of activated IP3 receptors. Mutations in this gene are a cause of spastic paraplegia-18 (SPG18). Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Feb 2012]

Product images:



Circular map for RC223513