

Product datasheet for RC223525L2

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn

OriGene Technologies, Inc.

GRIA1 (NM_000827) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: GRIA1 (NM_000827) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: GRIA1

Synonyms: GluA1; GLUH1; GLUR1; GLURA; HBGR1

Mammalian Cell None

Selection:

Vector:pLenti-C-mGFP (PS100071)E. coli Selection:Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC223525).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





 $[\]ensuremath{^*}$ The last codon before the Stop codon of the ORF.

ACCN: NM_000827

ORF Size: 2718 bp



GRIA1 (NM_000827) Human Tagged Lenti ORF Clone - RC223525L2

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 000827.2</u>

 RefSeq Size:
 3242 bp

 RefSeq ORF:
 2721 bp

 Locus ID:
 2890

 UniProt ID:
 P42261

Cytogenetics: 5q33.2

Domains:

Protein Families: Druggable Genome, Ion Channels: Glutamate Receptors, Transmembrane

lig chan, ANF receptor

Protein Pathways: Amyotrophic lateral sclerosis (ALS), Long-term depression, Long-term potentiation,

Neuroactive ligand-receptor interaction

MW: 101.3 kDa

Gene Summary: Glutamate receptors are the predominant excitatory neurotransmitter receptors in the

mammalian brain and are activated in a variety of normal neurophysiologic processes. These

receptors are heteromeric protein complexes with multiple subunits, each possessing transmembrane regions, and all arranged to form a ligand-gated ion channel. The


classification of glutamate receptors is based on their activation by different pharmacologic

agonists. This gene belongs to a family of alpha-amino-3-hydroxy-5-methyl-4-isoxazole propionate (AMPA) receptors. Alternatively spliced transcript variants encoding different

isoforms have been found for this gene. [provided by RefSeq, Jul 2008]



Product images:



Circular map for RC223525L2