

Product datasheet for **RC223532L3V**

Semaphorin 3B (SEMA3B) (NM_004636) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Semaphorin 3B (SEMA3B) (NM_004636) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Semaphorin 3B
Synonyms:	LUCA-1; SemA; SEMA5; SEMAA; semaV
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004636
ORF Size:	2247 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223532).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004636.2
RefSeq Size:	2989 bp
RefSeq ORF:	2250 bp
Locus ID:	7869
UniProt ID:	Q13214
Cytogenetics:	3p21.31
Domains:	Sema, IG, PSI
Protein Families:	Secreted Protein, Transmembrane



[View online »](#)

Protein Pathways: Axon guidance

MW: 83.12 kDa

Gene Summary: The protein encoded by this gene belongs to the class-3 semaphorin/collapsin family, whose members function in growth cone guidance during neuronal development. This family member inhibits axonal extension and has been shown to act as a tumor suppressor by inducing apoptosis. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Feb 2014]