

## Product datasheet for **RC223538L3V**

### **LYPD3 (NM\_014400) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	LYPD3 (NM_014400) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LYPD3
Synonyms:	C4.4A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_014400
ORF Size:	1038 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223538).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_014400.1</a>
RefSeq Size:	1745 bp
RefSeq ORF:	1041 bp
Locus ID:	27076
UniProt ID:	<a href="#">O95274</a>
Cytogenetics:	19q13.31
Domains:	LU
Protein Families:	Druggable Genome


[View online »](#)

**MW:** 36 kDa

**Gene Summary:** Supports cell migration. May be involved in urothelial cell-matrix interactions. May be involved in tumor progression.[UniProtKB/Swiss-Prot Function]