

Product datasheet for **RC223552L2V**

GPR128 (ADGRG7) (NM_032787) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GPR128 (ADGRG7) (NM_032787) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GPR128
Synonyms:	GPR128
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_032787
ORF Size:	2391 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223552).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032787.1
RefSeq Size:	3070 bp
RefSeq ORF:	2394 bp
Locus ID:	84873
UniProt ID:	Q96K78
Cytogenetics:	3q12.2
Protein Families:	Druggable Genome, GPCR, Transmembrane
MW:	88.7 kDa



[View online »](#)

Gene Summary: Orphan receptor.[UniProtKB/Swiss-Prot Function]