

Product datasheet for **RC223561L4V**

MAGEA2 (NM_005361) Human Tagged ORF Clone Lentiviral Particle

Product data:


Product Type:	Lentiviral Particles
Product Name:	MAGEA2 (NM_005361) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MAGEA2
Synonyms:	CT1.2; MAGE2; MAGEA2A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_005361
ORF Size:	942 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223561).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005361.2
RefSeq Size:	2024 bp
RefSeq ORF:	945 bp
Locus ID:	4101
UniProt ID:	P43356
Cytogenetics:	Xq28
MW:	35.1 kDa



[View online »](#)

Gene Summary:

This gene is a member of the MAGEA gene family. The members of this family encode proteins with 50 to 80% sequence identity to each other. The promoters and first exons of the MAGEA genes show considerable variability, suggesting that the existence of this gene family enables the same function to be expressed under different transcriptional controls. The MAGEA genes are clustered at chromosomal location Xq28. They have been implicated in some hereditary disorders, such as dyskeratosis congenita. This gene has two identical copies at different loci. Alternatively spliced transcript variants encoding the same protein have been identified for this gene. [provided by RefSeq, Jul 2008]

Product images:

[RC223561L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC223561L4V particle to overexpress human MAGEA2-mGFP fusion protein.