

Product datasheet for **RC223586**

GGT1 (NM_005265) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GGT1 (NM_005265) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	GGT1
Synonyms:	CD224; D22S672; D22S732; GGT; GGT 1; GTG
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC223586 representing NM_005265
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAAGAAGAAGTTAGTGGTCTGGCCTGCTGGCCGTGGTCTGGTCTGGTATTGTGCGCCTCTGTC
 TCTGGCTGCCCTCAGCCTCCAAGGAACCTGACAACCATGTGTACACCAGGGCTGCCGTGGCCGGATGC
 CAAGCAGTGCTCGAAGATTGGGAGGGATGCACTGCGGGACGGTGGCTCTGCGGTGGATGCAGCCATTGCA
 GCCCTGTTGTGTGGGGCTCATGAATGCCACAGCATGGGCATCGGGGGTGGCCTTCTCCTCACCATCT
 ACAACAGCACCACAGAAAAGCTGAGGTCAACACGCCCGCAGGTGGCCCCAGGCTGGCCTTTGCCAC
 CATGTTCAACAGCTCGGAGCAGTCCCAGAAGGGGGGCTGTGGTGGCGGTGCCTGGGGAGATCCGAGGC
 TATGAGCTGGCACACCAGCGCATGGGCGGTGCCCTGGGCTCGCCTCTCCAGCCCAGCATCCAGCTGG
 CCCGCCAGGGCTTCCCCGTGGGCAAGGGCTTGGCGGCAGCCCTGAAAAACAAGCGGACCGTCATCGAGCA
 GCAGCCTGTCTTGTGTGAGGTGTCTGCCGGGATAGAAAGGTGCTTCGGGAGGGGAGAGACTGACCCTG
 CCGCAGCTGGCTGACACCTACGAGACGCTGGCCATCGAGGGTGCCAGGCCTTACAACGGCAGCCTCA
 CGGCCAGATTGTGAAGGACATCCAGGCGGCCGGGGCATTGTGACAGCTGAGGACCTGAACAACCTACCG
 TGCTGAGCTGATCGAGCACCCGCTGAACATCAGCCTGGGAGACGTGGTGTGTACATGCCAGTGCGCCG
 CTCAGCGGGCCCGTGTGGCCCTCATCCTCAACATCCTCAAAGGGTACAACCTTCTCCCGGAGAGCGTGG
 AGAGCCCCGAGCAGAAGGGCCTGACGTACCACCGCATCGTAGAGGCTTCCGGTTTGCTACGCCAAGAG
 GACCCTGCTTGGGGACCCCAAGTTTGGGATGTGACTGAGGTGGTCCGCAACATGACCTCCGAGTTCTTC
 GCTGCCAGCTCCGGGCCAGATCTCTGACGACACCACTACCCGATCTCTACTACAAGCCCGAGTTCT
 ACACGCCGATGACGGGGCAGTCTCACCTGTCTGCTGCGCAGAGGACGGCAGTGTGTGTCGCCAC
 CAGCACCATCAACCTCTACTTTGGCTCCAAGGTCCGCTCCCGGTGAGCGGGATCCTGTTCAATAATGAA
 ATGGACGACTTACGCTCTCCAGCATACCAACGAGTTTGGGGTACCCCCCTCACCTGCCAATTTATCC
 AGCCAGGGAAGCAGCCGCTCTCGTCCATGTGCCGACGATCATGGTGGGCCAGGACGGCCAGGTCCGGAT
 GGTGGTGGGAGCTGTGGGGCACACAGATCACCACGGCCACTGCACTGGCCATCATCTACAACCTCTGG
 TTCGGCTATGACGTGAAGCGGGCGTGGAGGAGCCCGGCTGCACAACAGCTTCTGCCAACGTACGCA
 CAGTGGAGAGAAACATTGACCAGGCAGTACTGCAGCCCTGGAGACCCGGCACCATCACACCAGATCGC
 GTCCACCTTATCGCTGTGGTGAAGCCATCGTCCGACGGCTGGTGGCTGGGCAGCTGCCTCGGACTCC
 AGGAAAGCGGGGAGCCTGCCGGCTAC

ACGCGTACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC223586 representing NM_005265
 Red=Cloning site Green=Tags(s)

MKKKLVLVLLAVVLVIVVGLCLWLPSASKEPDNHVYTRAAVAADAKQCISKIGRDALRDGGSAYDAAIA
 ALLCVGLMNAHSMGIGGLFLTIYNSTTRKAEVINAREVAPRLAFATMFNSSEQSQKGLSVAVPGEIRG
 YELAHQRHGRLPWARLFQPSIQLARQGFVPGKGLAAALENKRTVIEQQPVLCEVFCRDRKVLREGERLTL
 PQLADTYETLAIEGAQAFYNGSLTAQIVKDIQAAGGIVTAEDLNRYRAELIEHPLNISLGDVVL YMP
 SAPLSGPVLAALILNLIKGYNFSRESVESPEQKGLTYHRIVEAFRFAYAKRTLLGDPKFVDVTEVVRNMTSEFF
 AAQLRAQISDDTTHPI SYKPEFYTPDDGGTAHL SVVAEDGSAVSATSTINLYFGSKVRSVPVSGILFNNE
 MDDFSSPSITNEFGVPPSPANFIQPGKQPLSSMCPTIMVGQDQVQRMVGAAGGTQITTALAI IYNLW
 FGYDVKRAVEEPRLNHQLLPNVTTVERNIDQAVTAALETRHHHTQIASTFIAVVQAVIRTAGGWAAASDS
 RKGGEPEGY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk6218_c09.zip

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:


ACCN: NM_005265

ORF Size: 1707 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_005265.3](#)

RefSeq Size: 2431 bp

RefSeq ORF: 1709 bp

Locus ID: 2678

Cytogenetics: 22q11.23

Domains: G_glu_transpept


Protein Families: Protease, Transmembrane

Protein Pathways: Arachidonic acid metabolism, Cyanoamino acid metabolism, Glutathione metabolism, Metabolic pathways, Selenoamino acid metabolism, Taurine and hypotaurine metabolism


MW: 61.2 kDa

Gene Summary: The enzyme encoded by this gene is a type I gamma-glutamyltransferase that catalyzes the transfer of the glutamyl moiety of glutathione to a variety of amino acids and dipeptide acceptors. The enzyme is composed of a heavy chain and a light chain, which are derived from a single precursor protein. It is expressed in tissues involved in absorption and secretion and may contribute to the etiology of diabetes and other metabolic disorders. Multiple alternatively spliced variants have been identified. There are a number of related genes present on chromosomes 20 and 22, and putative pseudogenes for this gene on chromosomes 2, 13, and 22. [provided by RefSeq, Jan 2014]

Product images:



Circular map for RC223586



Western blot validation of overexpression lysate (Cat# [LY422309]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with [RC223689] using transfection reagent MegaTran 2.0 (Cat# [TT210002]).