

Product datasheet for RC223589

Marapsin2 (PRSS38) (NM_183062) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Marapsin2 (PRSS38) (NM_183062) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Marapsin2
Synonyms:	MPN2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC223589 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGCCCTGCTTCCGTCATGGGCCACTCGGGCCCTCGCCCTGGGCTTCTGCTGCTGCTCCTGG
TGGTGGCCCTCCCGGGTCGCAGCATTGGTCCACAGACAGCCAGAGAACCAGGGAATCTCCCTAACTGG
CAGCGTGGCTGTGGTCGGCCAGCATGGAGGGAAAAATCCTGGGCGCGTCCCTGCGCCCAGAGGAAG
TGGCCGTGGCAGGTCAGCGTCACTACGCAGGCCCTCCACGTCTGCGGCGGCTCCATCCTCAATGAGTACT
GGGTGCTGTCAGCTGCGCACTGCTTTCACAGGGACAAGAATATCAAATCTATGACATGTACGTAGGCCCT
CGTAAACCTCAGGGTGGACGGCAACCACACCCAGTGGTATGAGGTGAACAGGGTGATCCTGCACCCACA
TATGAGATGTACCACCCATCGGAGGTGACGTGGCCCTGGTGCAGCTGAAGACCCGATTGTGTTTTCTG
AGTCCGTGCTCCCGGTTTGCCTTCAACTCCAGAAGTGAACCTTACCAGTGCCAATTGCTGGGCTACGGG
ATGGGGACTAGTCTAAAACAAGGTGAGACCTCAGACGAGCTGCAGGAGGTGCAGCTCCCGCTGATCCTG
GAGCCCTGGTGCCACCTGCTCTACGGACACATGCTCTACATCATGCCGACATGCTGTGTGCTGGGGACA
TCCTGAATGCTAAGACCGTGTGTGAGGGCGACTCCGGGGCCCACTTGTCTGTGAATTCAACCGCAGCTG
GTTGCAGATTGGAATTGTGAGCTGGGGCCGAGGCTGCTCAACCCCTGTACCCTGGAGTGTATGCCAGT
GTTTCCTATTTCTCAAATGGATATGTATAACATAGAAATCACGCCCACTCCTGCTCAGCCAGCCCTG
CTCTCTCCAGCTCTGGGGCCCACTCTCAGCGTCTAATGGCCATGCTGGCTGGCTGGTCAAGTCTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC223589 protein sequence
Red=Cloning site Green=Tags(s)

MAAPASVMGPLGPSALGLLLLLLVVAPPRVAALVHRQPENQGISLTGSVACGRPSMEGKILGGVPAPERK
 WPWQVSVHYAGLHVCGGSILNEYWVLSAAHCFHRDKNIKIYDMYVGLVNLRVLDGNHTQWYEVNRVILHPT
 YEMYHPIGGDVALVQLKTRIVFSESVLPVCLATPEVNLTSANCWATGWGLVSKQGETSDELQEVQLPLIL
 EPWCHLLYGHMSYIMPDMLCAGDILNAKTVCEGDSGGPLVCEFNRSWLQIGIVSWGRGCSNPLYPGVYAS
 VSYF SKWICDNIEITPTPAQPAPALSPALGPTLSVLMAMLAGWSVL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6605_g06.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_183062

ORF Size: 978 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_183062.2](#), [NP_898885.1](#)

RefSeq Size: 1243 bp

RefSeq ORF: 981 bp

Locus ID: 339501


UniProt ID: [A1L453](#)

Cytogenetics: 1q42.13


Protein Families: Druggable Genome

MW: 35.4 kDa

Product images:



Circular map for RC223589



Western blot validation of overexpression lysate (Cat# [LY405225]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC223589 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).