

Product datasheet for **RC223589L3V**

Marapsin2 (PRSS38) (NM_183062) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Marapsin2 (PRSS38) (NM_183062) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Marapsin2
Synonyms:	MPN2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_183062
ORF Size:	978 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223589).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_183062.2 , NP_898885.1
RefSeq Size:	1243 bp
RefSeq ORF:	981 bp
Locus ID:	339501
UniProt ID:	A1L453
Cytogenetics:	1q42.13
Protein Families:	Druggable Genome
MW:	35.4 kDa



[View online »](#)