

Product datasheet for **RC223611**

GIDRP88 (R3HCC1L) (NM_014472) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GIDRP88 (R3HCC1L) (NM_014472) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	GIDRP88
Synonyms:	C10orf28; GIDRP86; GIDRP88; PSORT
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC223611 ORF sequence

Red=Cloning site Blue=ORF Green=Tags(s)

 TTTTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

 ATGCAGCAAGAATCAGAGAGATGCAGAGTTAGAGCCAGAAGGCCTGACATGGCACTTTATGTACCTAAAG
CTCGTAGGGGTGCAGTACTCCTTAAGACAGGTGATGAAGAAGAAAGCTGTGGTTACCTAACTCTGTGGT
GAAAGAAAAGCAAAAAGAAAGTTCTCTCTCCAAAAAGAGTCTTTAAAGACAAACCGGAGGCTCGAAGA
CTAAATATCAATCCTGATAGAAAGGAGCATAATTGTAGAGAAGAAAAGAAATCTTCAACAAAATTAAGAA
TGGACACATGCCTTCAAAAAACAAATAGAGTTTGTCTAAGAGAGGAACCACTGAATCCAAAGAAGTATT
ATCCCAAGGACAACAGCAGGGAGCCCCAAATGCTGGGGTTATAACTAATGCACCTTTGCAGAGACATTTT
AAACCAAGAAGGTGGAGTGTGGGAAGTTGAACTACGGATGTGACAGGACATGAGAGGATACTTCTTT
CACAGGCTGTAGAAATCAGCGAGGCTCAAGTTCCAAGCAAACCATTCAAAAATGTGGAATTCTGTGA
CTTCAGTAGGCATGAACCTGATGGGAAGCATTTGAAGACAAAGATTTGGAAGGCAGAATTGAAACTGAT
ACCAAGGTTTTGGAGATACTATATGAGTTTCTAGAGTTTTAGTTCTGTGATGAAACCTGAGAATATGA
TTGTACCAATAAACTAAGCTCTGATTCTGAAATTGTACAACAAAGCATGCAACATCAGATGGAATATT
GAATCCCAGCAGCGGAGGCATCACCCTACTTCTGTTCTGGAAGTCCAGATGGTGTCTTTGATCAAACT
TGCGTAGATTTTGAAGTTGAGAGTGTAGGTGGTATAGCCAATAGTACAGGTTTCATCTTAGATCAAAAAG
ATACAGATTCCATTCTGCAACTATGGGTACATCTCTGTGACAGAGCACAATGACACTGTTAGTCC
AGTAATGATTAGAGAATGTGAGAAGAATGACAGCACTGCTGATGAGTTACATGTAAGCAGCAACCTCCT
GATACAGCTGTCTTGTCTATGAAACACATAGAGATAGTGGATTTAAGAATGTAGGTGACATTACCAATA
AAGCATGTATGATGGACACTACAGGTATGTCCTGTAGTGATCATGTAAGTGTGATAGCCCTTATGTAGT
TGCAGTTAGAATAGCTGATGAGACCTCTATTAATACACGAAGTTTCTCAAAGTTTGTAGGAATGAGTGCA
GATGCAACCCCTCTTCATGTAGCTAGAAGTGGGAATGACACTGAAGATTTGAGCAACCCCTCTGCTTGT
CAGATATTTATGGTGAGAGTATTTTCATCTCATTTTACAGAGTCAACAGGAAAGTTGATAGAGAGCTTGTC
AGATTGTGCTTCTCTTACCTATAAAAAAGATTGCTGGTAGTAATTATAACACTTTTTTGGACTCTGAA
CTCAGTATGTTAAATGGGACAAAAGTTCTTTCAGACAGTGCCGTGGGCATTGACCTGGGTAGTACTGGTG
ATACAACAGAAGCATTGCACGAATAAGAACTGCCGAAGAGTTCAAACAGAAGCAAGATGCCTCAGG
GAGTATAGAATTTGGTGTATCTTTCTGATAGGGAATCATCATCTATGGAACATCCATCGAACCAAAA
GCAACTGAACTTCTCACACAGAGGGAATTACTGCCATTGAGGAGAGCTGGGAGTCTATGTTTAACGATG
ATGGTGACTGCCTGGATCCACGTCTTCTACAAGAGTTATCAGGGAATACCAAGAGCAGAGAGAGCATCCA
GGAACCTAGATCTGATTACTACAATCATGAAGTTCTGATATTGACCTCAGTGATTGTGAATCCCACAT
GTCATTGAAATTTATGACTTTCCCAAGAATTTCTGACTGAAGACCTTCTACGGGTTTTCTGCAGTTATC
AAAAGAAAGGATTTGATATTAATGGGTGGATGATACACATGCCCTAGGAGTATTCTCCAGTCCAATTAC
AGCTCGTGATGCGTTGGGTATTAACACACCATGGTGAAGATTCGTCCCTTGTACAGGCCACAAGAGCA
GCCAAGGCCAAAGCTAGAGCTTATGCTGAGTTCTCCAGCCAGCAAAGGAGCGTCTGAGACTTCAGCAG
CCCTAGCCAGAAGTTAGTCATCAGTGCCCTTGGGGTTCGAAGTAAGCAGAGCAAAACCGAACGAGAAGC
AGAGCTCAAGAACTGCAAGAAGCCAGAGAGAGAAAGCGGTTGGAAGCCAAGCAACGGGAAGACATCTGG
GAAGGCAGAGACCAGTCTACAGTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC223611 protein sequence
 Red=Cloning site Green=Tags(s)


MQQESERCRVRARRPDMALYVPKARRGAVLLKTGDEEESCGSPNSVVKEKQKESLSQKEVFKDKPEAR
 LNINPDRKEHNCREEKKSSTKL RMDTCLQKTNRVCSKRGTTESKEVLSQGQQGAPNAGVITNAPLQRHF
 KPKKVECLEVETTDVTGHERILLSQACLEISEAQVPSKPFQNVFCDFSRHEPDGEAFEDKDLEGRIETD
 TKVLEILYEFPRVFSSVMKPENMIVPIKLSSDSEIVQQSMQTS DGILNPSSGGITTTSPVGPDPGVFDQT
 CVD FEVESVGGIANSTGFILDQKDTDSIPATMGHISLSESTNDTVSPVMIRECEKNDSTADELHVKEPP
 DTA VLAHETHRDSGFKNVGDITNKACMMDTTGMSCDHVTVDSPYVAVRIADETSINTRSF SKFVGMSA
 DATPLHVARSGNDTEDFS NPSACSDIYGESISSHFTTESTGKLIESLSDCASSLP IKKIAGSNYNTFLDSE
 LSMLNGTKVLSDSAVGIDL GSTGDTTEALHELRTAEFEKTEEQDASGSIEFGVSFPDRESSMETSI EPK
 ATETSHTEGITAEESWESMFNDGDCLDPRLQLSGNTKSRESIQEPRSDYYNHEVPDIDLSDCEFPH
 VIEIYDFPQEFRTEDLLRVFC SYQKGF DIKWDDTHALGVFSSPITARDALGIKHTMVKIRPLSQATRA
 AKAKARAYAEFLQPAKERPETS AALARRLVISALGVRSKQSKTEREAELKKLQEARERKRLEAKQREDIW
 EGRDQSTV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6554_b08.zip

Restriction Sites: Sgfi-MluI

Cloning Scheme:



ACCN: NM_014472

ORF Size: 2334 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_014472.5

RefSeq Size: 3386 bp


RefSeq ORF: 2337 bp

Locus ID: 27291


Cytogenetics: 10q24.2

MW: 86.3 kDa

Product images:



Circular map for RC223611



Western blot validation of overexpression lysate (Cat# [LY415254]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC223611 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).