

## Product datasheet for **RC223613L4V**

### Presenilin 2 (PSEN2) (NM\_000447) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Presenilin 2 (PSEN2) (NM_000447) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PSEN2
Synonyms:	AD3L; AD4; CMD1V; PS2; STM2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_000447
ORF Size:	1344 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223613).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000447.1</a>
RefSeq Size:	2236 bp
RefSeq ORF:	1347 bp
Locus ID:	5664
UniProt ID:	<a href="#">P49810</a>
Cytogenetics:	1q42.13
Domains:	Presenilin, PSN
Protein Families:	Druggable Genome, Protease, Transmembrane



[View online »](#)

**Protein Pathways:** Alzheimer's disease, Notch signaling pathway

**MW:** 50 kDa

**Gene Summary:** Alzheimer's disease (AD) patients with an inherited form of the disease carry mutations in the presenilin proteins (PSEN1 or PSEN2) or the amyloid precursor protein (APP). These disease-linked mutations result in increased production of the longer form of amyloid-beta (main component of amyloid deposits found in AD brains). Presenilins are postulated to regulate APP processing through their effects on gamma-secretase, an enzyme that cleaves APP. Also, it is thought that the presenilins are involved in the cleavage of the Notch receptor such that, they either directly regulate gamma-secretase activity, or themselves act as protease enzymes. Two alternatively spliced transcript variants encoding different isoforms of PSEN2 have been identified. [provided by RefSeq, Jul 2008]