

## Product datasheet for **RC223633L3V**

### SPG7 (NM\_199367) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SPG7 (NM_199367) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SPG7
Synonyms:	CAR; CMAR; PGN; SPG5C
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_199367
ORF Size:	1467 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223633).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_199367.1</a> , <a href="#">NP_955399.1</a>
RefSeq Size:	2328 bp
RefSeq ORF:	1470 bp
Locus ID:	6687
UniProt ID:	<a href="#">Q9UQ90</a>
Cytogenetics:	16q24.3
Protein Families:	Protease, Transmembrane
MW:	53.9 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a mitochondrial metalloprotease protein that is a member of the AAA family. Members of this protein family share an ATPase domain and have roles in diverse cellular processes including membrane trafficking, intracellular motility, organelle biogenesis, protein folding, and proteolysis. Mutations in this gene cause autosomal recessive spastic paraplegia 7. Two transcript variants encoding distinct isoforms have been identified. [provided by RefSeq, Mar 2014]