

Product datasheet for **RC223645A1V**

Human Ribosomal protein S11 (RPS11) (NM_001015) AAV Particle

Product data:

Product Type: AAV Particles
Tag: Myc-DDK
Symbol: Ribosomal protein S11
Synonyms: S11
Mammalian Cell None
Selection:
Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC223645 representing NM_001015
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**C

ATGGCGGACATTGAGCTGAGCGTGCCTACCAAAGCAGCCGACCATCTTTCAAACAAGAAGAGGGTCC
TGCTGGGAGAACTGGCAAGGAGAAGCTCCCGCGTACTACAAGAACATCGGTCTGGGCTTCAAGACACC
CAAGGAGGCTATTGAGGGCACCTACATTGACAAGAAATGCCCTTCACTGGTAATGTGTCCATTCGAGGG
CGGATCCTCTCTGGCGTGGTGACCAAGATGAAGATGCAGAGGACCATTTGTCATCCGCCGAGACTATCTGC
ACTACATCCGAAGTACAACCGCTTCGAGAAGCGCCACAAGAACATGTCTGTACACCTGTCCCCCTGCTT
CAGGGACGTCCAGATCGGTGACATCGTCACAGTGGCGAGTGCCGGCCTCTGAGCAAGACAGTGCGCTTC
AACGTGCTCAAGGTCACCAAGGCTGCCGGCACCAAGAAGCAGTTCCAGAAGTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC223645 representing NM_001015
Red=Cloning site Green=Tags(s)

MADIQTERAYQKQPTIFQNKRVLLGETGKEKLPRYYKNIGLGFKTPKEAIEGTYIDKKCPFTGNVSIRG
RILSGVVTMMQMRTIVIRRDYLHYIRKYNRFKRNHKNMSVHLSPCFRDVQIGDIVTVGECRPLSKTVRF
NVLKVTKAAGTKKQFQKF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



ACCN:	NM_001015
ORF Size:	474 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001015.3</u>
RefSeq Size:	641 bp
RefSeq ORF:	477 bp
Locus ID:	6205
UniProt ID:	<u>P62280</u>
Cytogenetics:	19q13.33
MW:	18.3 kDa