

## Product datasheet for **RC223657L2V**

### HS2ST1 (NM\_012262) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	HS2ST1 (NM_012262) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HS2ST1
Synonyms:	dj604K5.2; NFSRA
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_012262
ORF Size:	1068 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223657).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012262.2</a>
RefSeq Size:	6708 bp
RefSeq ORF:	1071 bp
Locus ID:	9653
UniProt ID:	<a href="#">Q7LGA3</a>
Cytogenetics:	1p22.3
Domains:	HS2ST
Protein Pathways:	Heparan sulfate biosynthesis



[View online »](#)

**MW:** 41.7 kDa

**Gene Summary:** Heparan sulfate biosynthetic enzymes are key components in generating a myriad of distinct heparan sulfate fine structures that carry out multiple biologic activities. This gene encodes a member of the heparan sulfate biosynthetic enzyme family that transfers sulfate to the 2 position of the iduronic acid residue of heparan sulfate. The disruption of this gene resulted in no kidney formation in knockout embryonic mice, indicating that the absence of this enzyme may interfere with the signaling required for kidney formation. Two alternatively spliced transcript variants that encode different proteins have been found for this gene. [provided by RefSeq, Aug 2008]