

## Product datasheet for **RC223702L3V**

### Angiotensin Converting Enzyme 1 (ACE) (NM\_000789) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Angiotensin Converting Enzyme 1 (ACE) (NM_000789) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ACE
Synonyms:	ACE1; CD143; DCP; DCP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_000789
ORF Size:	3918 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223702).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000789.3</a>
RefSeq Size:	4969 bp
RefSeq ORF:	3921 bp
Locus ID:	1636
UniProt ID:	<a href="#">P12821</a>
Cytogenetics:	17q23.3
Domains:	Peptidase_M2



[View online »](#)

<b>Protein Families:</b>	Druggable Genome, ES Cell Differentiation/IPS, Protease, Secreted Protein, Transmembrane
<b>Protein Pathways:</b>	Hypertrophic cardiomyopathy (HCM), Renin-angiotensin system
<b>MW:</b>	149.7 kDa
<b>Gene Summary:</b>	<p>This gene encodes an enzyme involved in blood pressure regulation and electrolyte balance. It catalyzes the conversion of angiotensin I into a physiologically active peptide angiotensin II. Angiotensin II is a potent vasopressor and aldosterone-stimulating peptide that controls blood pressure and fluid-electrolyte balance. This angiotensin converting enzyme (ACE) also inactivates the vasodilator protein, bradykinin. Accordingly, the encoded enzyme increases blood pressure and is a drug target of ACE inhibitors, which are often prescribed to reduce blood pressure. This enzyme additionally plays a role in fertility through its ability to cleave and release GPI-anchored membrane proteins in spermatozoa. Many studies have associated the presence or absence of a 287 bp Alu repeat element in this gene with the levels of circulating enzyme. This polymorphism, as well as mutations in this gene, have been implicated in a wide variety of diseases including cardiovascular pathophysiologies, psoriasis, renal disease, stroke, and Alzheimer's disease. Regulation of the homologous ACE2 gene may be involved in progression of disease caused by several human coronaviruses, including SARS-CoV and SARS-CoV-2. Alternative splicing results in multiple transcript variants encoding both somatic (sACE) and male-specific testicular (tACE) isoforms. [provided by RefSeq, Sep 2020]</p>