

Product datasheet for **RC223711A1V**

Human TA3 (TAAR9) (NM_175057) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TA3 (TAAR9) (NM_175057) AAV Particle
Tag:	Myc-DDK
Symbol:	TA3
Synonyms:	TA3; TAR3; TAR9; TRAR3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC223711 representing NM_175057 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGTGAACAATTTCTCCCAAGCTGAGGCTGTGGAGCTGTGTTACAAGACGTGAACGAATCCTGCATTA
AAACTCCTTACTCGCCAGGTCCTCGATCTATCCTCTACGCCGTCCTTGGTTTTGGGGCTGTGCTGGCAGC
GTTTGAAACTTACTGGTCATGATTGCTATCCTTCACTTCAAACAACGCACACACCTACAACTTTCTG
ATTGCGTCGCTGGCCTGTGCTGACTTCTTGGTGGGAGTCACTGTGATGCCCTTCAGCACAGTGAGGCTG
TGGAGAGCTGTTGGTACTTTGGGACAGTTACTGTAAATCCATACATGTTTTGACACATCCTTCTGTTT
TGCTTCTTTATTTCAATTTATGCTGTATCTCTGTTGATAGATACATTGCTGTTACTGATCCTCTGACCTAT
CCAACCAAGTTTACTGTGTCAGTTTCAGGGATATGCATTGTTCTTTCTGGTTCTTTCTGTCACATACA
GCTTTTCGATCTTTTACACGGGAGCCAACGAAGAAGGAATTGAGGAATTAGTAGTTGCTCTAACCTGTGT
AGGAGGCTGCCAGGCTCCACTGAATCAAACCTGGGTCCTACTTTGTTTCTTCTATTCTTTATACCCAAT
GTCGCCATGGTGTATATACAGTAAGATATTTTGGTGGCCAAGCATCAGGCTAGGAAGATAGAAAGTA
CAGCCAGCCAAGCTCAGTCCTCCTCAGAGAGTTACAAGGAAAGAGTAGCAAAAAGAGAGAGAAAGGCTGC
CAAAACCTTGGAATTGCTATGGCAGCATTTCTGTCTCTTGGCTACCATACCTCGTTGATGCAGTGATT
GATGCTTATATGAATTTATACTCCTCTTATGTTTATGAGATTTTAGTTTGGTGTGTTTATTATAATT
CAGCTATGAACCCCTTGATTTATGCTTTCTTTACCAATGGTTTGGGAAGGCAATAAACTTATTGTAAG
CGGCAAGGTCTTAAGGACTGATTCGTCAACAATAATTTATTTCTGAAGAAGTAGAGACAGAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC223711 representing NM_175057 Red=Cloning site Green=Tags(s) MVNNFSQAEAVELCYKNVNESCIKTPYSPGPRSILYAVLGFGAVLAAFGNLLVMIAILHFKQLHTPTNFI IASLACADFLVGVTVMFPSTVRSVESCWYFGDSYCKFHTCFDTSFCFASLFHLCCISVDRIYAVTDPLTY PTKFTVSVSGICIVLSWFFSVTYSFSIFYTGANEEGIEELVVALTCVGGCQAPLNQNWLLCFLFFIPN VAMVFIYSKIFLVAKHQARKIESTASQAQSSSESYKERVAKRERKAAKTLGIAMAAFLVSWLPYLVDAVI DAYMNFITPPYVYEILVWCYYNSAMNPLIYAFFYQWFGKAIKLIVSGKVLRTDSSTTNLFSEEVETD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_175057
ORF Size:	1044 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_175057.3, NP_778227.3</u>
RefSeq Size:	1049 bp
RefSeq ORF:	1047 bp
Locus ID:	134860
UniProt ID:	<u>Q96RI9</u>
Cytogenetics:	6q23.2
MW:	38.8 kDa