

Product datasheet for **RC223740L3V**

USP44 (NM_001042403) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | USP44 (NM_001042403) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | USP44 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001042403 |
| ORF Size: | 2136 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC223740). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001042403.1 |
| RefSeq Size: | 3334 bp |
| RefSeq ORF: | 2139 bp |
| Locus ID: | 84101 |
| UniProt ID: | Q9H0E7 |
| Cytogenetics: | 12q22 |
| Protein Families: | Druggable Genome, Protease |
| MW: | 81 kDa |



[View online »](#)

Gene Summary:

The protein encoded by this gene is a protease that functions as a deubiquitinating enzyme. The encoded protein is thought to help regulate the spindle assembly checkpoint by preventing early anaphase onset. This protein specifically deubiquitinates CDC20, which stabilizes the anaphase promoting complex/cyclosome. [provided by RefSeq, Dec 2016]