

Product datasheet for **RC223751L3V**

UNC45B (NM_173167) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	UNC45B
Synonyms:	CMYA4; CTRCT43; MFM11; SMUNC45; UNC-45B; UNC45
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_173167
ORF Size:	2793 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC223751).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_173167.2
RefSeq Size:	5694 bp
RefSeq ORF:	2796 bp
Locus ID:	146862
UniProt ID:	Q8IWX7
Cytogenetics:	17q12



MW: 103.6 kDa

Gene Summary:

This gene encodes a co-chaperone required for folding and accumulation of type II myosins. The protein consists of three tetratricopeptide repeat motifs at the N-terminus that form a complex with heat shock protein 90, a central region of unknown function that is conserved in all Unc-45 proteins, and a C-terminal Unc-45/Cro1/She4 domain. The protein is expressed at high levels in striated muscle, where its muscle myosin chaperone activity is dependent on heat shock protein 90 acting as a co-chaperone. A missense mutation in this gene has been associated with cataract development. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2015]