

Product datasheet for RC223781A1V

Human H4C14 (NM_003548) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human H4C14 (NM_003548) AAV Particle
Tag:	Myc-DDK
Symbol:	H4C14
Synonyms:	FO108; H4; H4-16; H4/n; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4C11; H4C12; H4C13; H4C15; H4F2; H4FN; HIST2H4; HIST2H4A
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC223781 representing NM_003548 Red=Cloning site Blue=ORF Green=Tags(s) TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCCGCGATCGCC ATGTCCGGCAGAGGAAAGGGCGGAAAAGGCTTAGGCAAAGGGGGCGCTAAGCGCCACCGCAAGGTCTTGA GAGACAACATTCAGGGCATCACCAAGCCTGCCATTCGGCGTCTAGCTCGGCGTGGCGGCGTTAAGCGGAT CTCTGGCCTCATTTACGAGGAGACCCGCGGTGTGCTGAAGGTGTTCTCGGAGAATGTGATTCTGGGACGCA GTCACCTACACCGAGCACGCCAAGCGCAAGACCGTCACAGCCATGGATGTGGTGTACGCGCTCAAGCGCC AGGGGCGCACCCCTGTACGGCTTCGGAGGC ACGCGTACGCGGCGGCTCGAGCAGAAACTCATCTCTGAAGAGGATCTGGCAGCAATGATATCCTGGATT ACAAGGATGACGACGATAAGGTTTAA >RC223781 representing NM_003548 Red=Cloning site Green=Tags(s) MSGRGKGGKGLGKGGAKRHRKVLRDNIQGITKPAIRRLARRGGVKRISGLIYEETRGVLKVFLENVIRDA VTYTEHAKRKTVTAMDVVYALKRQGRPLYGFGG* TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Protein Sequence:	
Species:	Human
Serotype:	AAV-2
ACCN:	NM_003548



[View online »](#)

ORF Size:	309 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_003548.2</u> , <u>NP_003539.1</u>
RefSeq Size:	396 bp
RefSeq ORF:	312 bp
Locus ID:	8370
UniProt ID:	<u>P62805</u>
Cytogenetics:	1q21.2
MW:	11.2 kDa