

Product datasheet for **RC223804A1V**

Human Clorf84 (NM_001012961) AAV Particle

Product data:

Product Type: AAV Particles

Tag: Myc-DDK

Symbol: Clorf84

Synonyms: FLJ34502; RP11-506B15.1; SZT2A

Mammalian Cell None

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC223804 representing NM_001012961
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCTCGGAGCGCCCGGAGCCGGAGGTGGAAGAAGCTGGGCAGGTGTTCTGTAAATGAAAAAGGATT
ATCGAATCTCCGAAATGTTGCCTGGCTTGGTTCCTCAGTCATCTGCACAACTGTGCAGGCCACACC
CCAGGAGATGCTGTCTGAACAGGAATTGGAAGTCCTCAGTGTCTGCCCTGGGTGGCAGCCAGATGAA
CCAGTGGTCCCAAGGCCATTCTCTGGTACCTTCACCCGGGTACCTTCCTGGCTTGGCAGTATCGGT
TTGTCATTGAGTTGGACCTTAGCCCATCTACTGGCATTGTGGATGATCCACAGGGGAGATCTTGTGGA
TGAAGTTTTCCATGCCCTGTCCCGTGTCTAGGCGGGCTGCTTCGGCCCTTCCGAGTGCCTGGATCTTGC
ATCGACTTCCAGCCTGAGATCTATGTAACATCCAGGCCTACTCCTCCATCATTGGACTGCAGTCCCACC
AGAGACAGGGTTTCACCATGTTGGCCAGGCTGCCTCGAACTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTAA

Protein Sequence: >RC223804 representing NM_001012961
Red=Cloning site Green=Tags(s)

MASERPEPEVEEAGQVFLMKDYRISRNVLAWFLSHLHQTVPQATPQEMLSEQEVLVSVLPPGWQPD
PVVPRPFLVPSTRVTFWQYRFVIELDSPSTGIVDDSTGEILFDEVFHALSRCLGGLLRPFRVPGSC
IDFQPEIYVTIQAYSSIIGLQSHQRQGFMTLARLPSNF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2



ACCN:	NM_001012961
ORF Size:	534 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001012961.1</u> , <u>NP_001012979.1</u>
RefSeq Size:	1213 bp
RefSeq ORF:	536 bp
Locus ID:	149469
Cytogenetics:	1p34.2
MW:	20.1 kDa