

Product datasheet for **RC223814L4V**

MRC1L1 (NM_001009567) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MRC1L1 (NM_001009567) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MRC1L1
Synonyms:	bA541119.1; CLEC13DL
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001009567
ORF Size:	4368 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223814).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001009567.1 , NP_001009567.1
RefSeq Size:	5171 bp
RefSeq ORF:	4370 bp
Locus ID:	414308
Cytogenetics:	10p12.33
Protein Families:	Druggable Genome, Transmembrane
MW:	166 kDa



[View online »](#)

Gene Summary:

This gene encodes a type I membrane receptor protein highly related to the mannose receptor C type 1 protein (MRC1). The two genes are located near each other on chromosome 10 in a region thought to represent a segmental duplication. [provided by RefSeq]