

## Product datasheet for **RC223815L4V**

### CASC1 (CFAP94) (NM\_018272) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CASC1 (CFAP94) (NM_018272) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CFAP94
Synonyms:	CASC1; LAS1; PPP1R54
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018272
ORF Size:	2166 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223815).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018272.3</a>
RefSeq Size:	2523 bp
RefSeq ORF:	2169 bp
Locus ID:	55259
UniProt ID:	<a href="#">Q6TDU7</a>
Cytogenetics:	12p12.1
MW:	83.6 kDa



[View online »](#)

**Gene Summary:**

Via its association with the multisubunit axonemal dynein complex, is potentially involved in the regulation of cilia function. May also act as a cell cycle regulator.[UniProtKB/Swiss-Prot Function]