

Product datasheet for RC223820L1

SIAH3 (NM_198849) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Tag: Myc-DDK

Symbol: SIAH3

Mammalian Cell None

Selection:

Vector: pLenti-C-Myc-DDK (PS100064)



E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC223820).

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN: NM_198849

ORF Size: 807 bp




OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_198849.2
RefSeq Size:	4005 bp
RefSeq ORF:	810 bp
Locus ID:	283514
UniProt ID:	Q8IW03
Cytogenetics:	13q14.13
MW:	30.5 kDa
Gene Summary:	Negative regulator of PRKN translocation to damaged mitochondria. Acts probably by destabilizing PINK1 protein, hence inhibiting PRKN targeting to dysfunctional depolarized mitochondria.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC223820L1



Double digestion of RC223820L1 using SgfI and MluI