

Product datasheet for **RC223829L3V**

FUT4 (NM_002033) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FUT4 (NM_002033) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FUT4
Synonyms:	CD15; ELFT; FCT3A; FUC-TIV; FUTIV; LeX; SSEA-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002033
ORF Size:	1590 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223829).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002033.2
RefSeq Size:	5937 bp
RefSeq ORF:	1593 bp
Locus ID:	2526
UniProt ID:	P22083
Cytogenetics:	11q21
Protein Families:	Transmembrane
Protein Pathways:	Glycosphingolipid biosynthesis - lacto and neolacto series, Metabolic pathways



[View online »](#)

MW: 58.9 kDa

Gene Summary: The product of this gene transfers fucose to N-acetyllactosamine polysaccharides to generate fucosylated carbohydrate structures. It catalyzes the synthesis of the non-sialylated antigen, Lewis x (CD15). [provided by RefSeq, Jan 2009]