

Product datasheet for **RC223857**

Inositol Hexakisphosphate Kinase 2 (IP6K2) (NM_001005911) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Inositol Hexakisphosphate Kinase 2 (IP6K2) (NM_001005911) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: IP6K2
Synonyms: IHPK2; InsP6K2; PIUS
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC223857 representing NM_001005911
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAGCCAGCCTTCAGGGCCATGGATGTGGAGCCCCGCGCCAAAGGCGTCCTTCTGGAGCCCTTTGTCC
ACCAGGTCGGGGGCACTCATGCGTGCTCCGCTTCAATGAGACAACCCTGTGCAAGCCCTGGTCCCAAG
GGAACATCAGTTCTACGAGACCCTCCCTGCTGAGATGCGCAAATTCACTCCCCAGTACAAAGGACAAAGC
CAAAGGCCCTTGTAGCTGGCCATCCCTGCCCAATTTTTCCCTGGTCTTTCCCTGTGCCACAGG
GAAGTGTGGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC223857 representing NM_001005911
Red=Cloning site Green=Tags(s)

MSPAFRAMDVEPRAKGVLLPEPFVHQVGGHSCVLRFNETTLCCKPLVPREHQFYETLPAEMRKFTPQYKQGS
QRPLVSWPSLPHFFPWSFPLWPQGSVA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001005911

ORF Size: 291 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001005911.3
RefSeq Size:	1254 bp
RefSeq ORF:	294 bp
Locus ID:	51447
UniProt ID:	Q9UHH9
Cytogenetics:	3p21.31
Protein Families:	Druggable Genome
MW:	11.2 kDa
Gene Summary:	This gene encodes a protein that belongs to the inositol phosphokinase (IPK) family. This protein is likely responsible for the conversion of inositol hexakisphosphate (InsP6) to diphosphoinositol pentakisphosphate (InsP7/PP-InsP5). It may also convert 1,3,4,5,6-pentakisphosphate (InsP5) to PP-InsP4 and affect the growth suppressive and apoptotic activities of interferon-beta in some ovarian cancers. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]