

## Product datasheet for **RC223869L3V**

### CREB5 (NM\_004904) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CREB5 (NM_004904) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CREB5
Synonyms:	CRE-BPA; CREB-5; CREBPA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004904
ORF Size:	1503 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223869).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004904.2</a>
RefSeq Size:	8213 bp
RefSeq ORF:	1506 bp
Locus ID:	9586
UniProt ID:	<a href="#">Q02930</a>
Cytogenetics:	7p15.1
Domains:	BRLZ, zf-C2H2
Protein Families:	Transcription Factors



[View online »](#)

**Protein Pathways:** Huntington's disease, Prostate cancer

**MW:** 55.9 kDa

**Gene Summary:** The product of this gene belongs to the CRE (cAMP response element)-binding protein family. Members of this family contain zinc-finger and bZIP DNA-binding domains. The encoded protein specifically binds to CRE as a homodimer or a heterodimer with c-Jun or CRE-BP1, and functions as a CRE-dependent trans-activator. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]