

## Product datasheet for **RC223875L3V**

### GPR17 (NM\_005291) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	GPR17 (NM_005291) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GPR17
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_005291
ORF Size:	1017 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223875).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005291.1</a>
RefSeq Size:	2070 bp
RefSeq ORF:	1104 bp
Locus ID:	2840
UniProt ID:	<a href="#">Q13304</a>
Cytogenetics:	2q14.3
Protein Families:	Druggable Genome, GPCR, Transmembrane
MW:	37.86 kDa



[View online »](#)

**Gene Summary:**

Dual specificity receptor for uracil nucleotides and cysteinyl leukotrienes (CysLTs). Signals through G(i) and inhibition of adenylyl cyclase. May mediate brain damage by nucleotides and CysLTs following ischemia.[UniProtKB/Swiss-Prot Function]