

Product datasheet for **RC223883L4V**

CGGBP1 (NM_001008390) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | CGGBP1 (NM_001008390) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | CGGBP1 |
| Synonyms: | CGGBP; p20-CGGBP |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001008390 |
| ORF Size: | 501 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC223883). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001008390.1 |
| RefSeq Size: | 4608 bp |
| RefSeq ORF: | 504 bp |
| Locus ID: | 8545 |
| UniProt ID: | Q9UFW8 |
| Cytogenetics: | 3p11.1 |
| MW: | 18.8 kDa |


[View online »](#)

Gene Summary:

This gene encodes a CGG repeat-binding protein that primarily localizes to the nucleus. CGG trinucleotide repeats are implicated in many disorders as they often act as transcription- and translation-regulatory elements, can produce hairpin structures which cause DNA replication errors, and form regions prone to chromosomal breakage. CGG repeats are also targets for CpG methylation. In addition to its ability to bind CGG repeats and regulate transcription, this gene is believed to play a role in DNA damage repair and telomere protection. In vitro studies indicate this protein does not bind to methylated CpG sequences. [provided by RefSeq, Jul 2017]