

## Product datasheet for **RC223890L4V**

### USP28 (NM\_020886) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	USP28 (NM_020886) Human Tagged ORF Clone Lentiviral Particle
Symbol:	USP28
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020886
ORF Size:	3234 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223890).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_020886.2</a>
RefSeq Size:	4642 bp
RefSeq ORF:	3234 bp
Locus ID:	57646
UniProt ID:	<a href="#">Q96RU2</a>
Cytogenetics:	11q23.2
Domains:	UCH
Protein Families:	Protease
MW:	122.3 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a deubiquitinase involved in the DNA damage pathway and DNA damage-induced apoptosis. Overexpression of this gene is seen in several cancers. [provided by RefSeq, Oct 2016]