

Product datasheet for **RC223910L3V**

SMIM15 (NM_001048249) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	SMIM15
Synonyms:	C5orf43
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001048249
ORF Size:	222 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC223910).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001048249.1
RefSeq Size:	2894 bp
RefSeq ORF:	225 bp
Locus ID:	643155
UniProt ID:	Q7Z3B0
Cytogenetics:	5q12.1



Protein Families: Transmembrane

MW: 8.4 kDa