

Product datasheet for **RC223911L3V**

AHDC1 (NM_001029882) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	AHDC1
Synonyms:	MRD25; XIGIS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001029882
ORF Size:	4809 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC223911).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001029882.1 , NP_001025053.1
RefSeq Size:	5714 bp
RefSeq ORF:	4812 bp
Locus ID:	27245
UniProt ID:	Q5TGY3
Cytogenetics:	1p36.11-p35.3



MW: 168.2 kDa

Gene Summary: This gene encodes a protein containing two AT-hooks, which likely function in DNA binding. Mutations in this gene were found in individuals with Xia-Gibbs syndrome. [provided by RefSeq, Jun 2014]