

Product datasheet for **RC223914L3V**

FAM26E (CALHM5) (NM_153711) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	CALHM5
Synonyms:	C6orf188; dJ493F7.3; FAM26E
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_153711
ORF Size:	927 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC223914).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_153711.1
RefSeq Size:	1777 bp
RefSeq ORF:	930 bp
Locus ID:	254228
UniProt ID:	Q8N5C1
Cytogenetics:	6q22.1



Protein Families: Transmembrane

MW: 35 kDa

Gene Summary: Pore-forming subunit of a voltage-gated ion channel.[UniProtKB/Swiss-Prot Function]